**MIRASOL®**

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**A Novel Approach to Blood Safety**

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## References

1. Brookmeyer R, Goedert JJ. Censoring in an epidemic with an application to hemophilia-associated AIDS. *Biometrics*. 1989;45:325-335.
2. Janatpour KA, Holland PV. A brief history of blood transfusion. In: Hillier CD, Silberstein LE, Ness PM, et al., eds. *Blood Banking and Transfusion Medicine: Basic Principles and Practice*. 2<sup>nd</sup> ed. Philadelphia, PA: Churchill Livingstone Elsevier; 2007:3-11.
3. Havig T. Where will pathogen inactivation have the greatest impact? *ISBT Science Series*. 2007;2:25-29.
4. Chamberland ME. Emerging infectious agents: do they pose a risk to the safety of transfused blood and blood products? *Clin Infect Dis*. 2002;34:797-805.
5. Soldan K, Barbara J. The risks of infection transmission by blood transfusions in England. *J Clin Pathol*. 1999;52:405-408.
6. Blajchman MA, Beckers EA, Dickmeiss E, et al. Bacterial detection of platelets: current problems and possible resolutions. *Transfus Med Rev*. 2005;19:259-272.
7. Chiu EKW, Yuen KY, Lie AKW, et al. A prospective study of symptomatic bacteremia following platelet transfusion and of its management. *Transfusion*. 1994;34:950-954.
8. Sullivan MT, McCullough J, Schreiber GB, Wallace EL. Blood collection and transfusion in the United States in 1997. *Transfusion*. 2002;42:1253-1260.
9. Andreu G, Morel P, Forestier F, et al. Haemovigilance network in France: organization and analysis of immediate transfusion incident reports from 1994 to 1998. *Transfusion*. 2002;42:1356-1364.
10. The Trial to Reduce Alloimmunization to Platelets Study Group. Leukocyte reduction and ultraviolet B irradiation of platelets to prevent alloimmunization and refractoriness to platelet transfusions. *N Engl J Med*. 1997;337:1861-1869.
11. Andreu G, Dewailly J. Prevention of HLA alloimmunization by using leukocyte-depleted components. *Curr Stud Hematol Blood Transfus*. 1994;60:29-40.
12. Anderson NA, Gray S, Copplestone JA, et al. A prospective randomized study of three types of platelet concentrates in patients with haematological malignancy: corrected platelet increments and frequency of non haemolytic febrile transfusion reactions. *Transfus Med*. 1997;7:33-39.
13. Enright H, Davis K, Gernsheimer T, et al. Factors influencing moderate to severe reactions to PLT transfusions: experience of the TRAP multicenter clinical trial. *Transfusion*. 2003;43:1545-1552.
14. Thaler M, Shamiss A, Orgad S, et al. The role of the blood from HLA-homozygous donors in fatal transfusion associated graft-versus-host disease after open heart surgery. *N Engl J Med*. 1989;321:25-28.
15. Wagner FF, Flegel WA. Transfusion-associated graft-versus host disease: risk due to homozygous HLA haplotypes. *Transfusion*. 1995;35:284-291.
16. Seghatchian J, de Sousa G. Pathogen-reduction systems for blood components: the current position and future trends. *Transfus Apher Sci*. 2006;35:189-196.
17. Pelletier JPR, Transue S, Snyder EL. Pathogen inactivation techniques. *Best Pract Res Clin Haematol*. 2006;19:205-242.
18. Moeremans K, Warie H, Annemans L. Assessment of the economic value of the INTERCEPT blood system in Belgium. *Transfus Med*. 2006;16:17-30.
19. Weusten JJ, Van Drimelen HA, Lelie PN. Mathematic modeling of the risk of HBV, HCV and HIV transmission by window-phase donations not detected by NAT. *Transfusion*. 2002;42:537-548.
20. Alter HJ, Stramer SL, Dodd RY. Emerging infectious diseases that threaten the blood supply. *Semin Hematol*. 2007;44:32-41.
21. Watson R. Europe witnesses first local transmission of chikungunya fever in Italy. *BMJ*. 2007;335:532-533.
22. Eder AF, Kennedy JM, Dy BA, et al. American Red Cross Regional Blood Centers. Bacterial screening of apheresis platelets and the residual risk of septic transfusion reactions: the American Red Cross experience (2004-2006). *Transfusion*. 2007;47:1134-1142.
23. Akahoshi M, Takanashi M, Masuda M, et al. A case of transfusion-associated graft-versus-host disease not prevented by white cell-reduction filters. *Transfusion*. 1992;32:169-172.
24. Visconti MR, Pennington J, Garner SF, et al. Assessment of removal of human cytomegalovirus from blood components by leukocyte depletion filters using real-time quantitative PCR. *Blood*. 2004;103:1137-1139.
25. Pennington J, Taylor GP, Sutherland J, et al. Persistence of HTLV-I in blood components after leukocyte depletion. *Blood*. 2002;100:677-681.
26. Roberts P. Resistance of vaccinia virus to inactivation by solvent detergent treatment of blood products. *Biologicals*. 2000;28:29-32.

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## References

1. Kumar V, Lockertie O, Keil SD, et al. Riboflavin and UV-light based pathogen reduction: extent and consequence of DNA damage at the molecular level. *Photochem Photobiol*. 2004;80:15–21.
2. World Health Organization. Thiamine, riboflavin, niacin, vitamin B6, pantothenic acid, and biotin. In: *Vitamin and Mineral Requirements in Human Nutrition*. 2nd ed. 2004;Section 9.3.3:172.
3. Scientific Committee for Food (SCF) (2000). Tolerable upper intake level of vitamin B2. (Opinion expressed on 22 November 2000.)
4. Rivlin SR. Riboflavin metabolism. *N Engl J Med*. 1970;285:463–472.
5. Kale H, Hankumar P, Nair PM, Netravali MS. Assessment of the genotoxic potential of riboflavin and lumiflavin. *Mutat Res*. 1992;298:9–16.
6. United States Food and Drug Administration (FDA) Center for Food Safety and Applied Nutrition (CFSAN) Priority-based Assessment of Food Additives (PAFA). EAFUS: A Food Additive Database. Available at <http://vm.cfsan.fda.gov/~7Edms/efalus.html>. Accessed November 2007.
7. Warburg VO, Christian W. Über das neue Oxydationsstament. *Naturwissenschaften*. 1932;20:990–991.
8. Speck WT, Rosenkranz S, Rosenkranz HS. Further observations on the photo-oxidation of DNA in the presence of riboflavin. *Biochim Biophys Acta*. 1976;435:39–44.
9. Oster G, Bellin JS, Holmstrom B, et al. Photochemistry of riboflavin. *Experientia*. 1962;18:249–253.
10. Tsujita A, Okada Y, Uehara K. Photosensitized inactivation of ribonucleic acids in the presence of riboflavin. *Biochim Biophys Acta*. 1965;103:360–363.
11. Speck WT, Chen CC, Rosenkranz HS. In vitro studies of effects of light and riboflavin on DNA and HeLa cells. *Pediatr Res*. 1975;9:150–153.
12. Kuratomi K, Kobayashi Y. Studies on the interactions between DNA and flavins. *Biochim Biophys Acta*. 1977;476:207–217.
13. Korycka-Dahl M, Richardson T. Photodegradation of DNA with fluorescent light in the presence of riboflavin and photoprotection by flavin triplet-state quenchers. *Biochim Biophys Acta*. 1990;610:229–230.
14. Cairns WL, Metzler DE. Photochemical degradation of flavins. VI. A new photoproduct and its use in studying the photolytic mechanism. *J Am Chem Soc*. 1971;93:2772–2777.
15. Hardwick CC, Herivel TR, Hernandez SC, et al. Separation, identification and quantification of riboflavin and its photoproducts in blood products using high-performance liquid chromatography with fluorescence detection: a method to support pathogen reduction technology. *Photochem Photobiol*. 2004;80:609–615.
16. Sisson TR. Photodegradation of riboflavin in neonates. *Federation Proc*. 1987;46:1883–1885.
17. Olsen JH, Hertz H, Kjaer SK, et al. Childhood leukemia following phototherapy for neonatal hyperbilirubinemia (Denmark). *Cancer Causes Control*. 1996;7:411–414.
18. Goodrich RP, Edrich RA, Goodrich LL, et al. The antiviral and antibacterial properties of riboflavin and light: applications to blood safety and transfusion medicine. In: *Comprehensive Series in Photochemical & Photobiological Sciences*. 2006;6:83–113.
19. Martin CB, Willing E, Ruane P, et al. An action spectrum of the riboflavin-photosensitized inactivation of lambda phage. *Photochem Photobiol*. 2005;81:474–480.
20. Dardare H, Platz MS. Binding affinities of commonly employed sensitizers of viral inactivation. *Photochem Photobiol*. 2002;75:561–564.
21. Joshi PD. Ultraviolet radiation-induced photodegradation and  $^{18}O$ ,  $O_2$  production by riboflavin, lumichrome and lumiflavin. *Indian J Biochem Biophys*. 1989;26:186–189.
22. Data on file. Navigant Biotechnologies LLC.
23. Goodrich RP, Edrich RA, Li J, et al. The Mirasol™ PRT system for pathogen reduction of platelets and plasma: an overview of current status and future trends. *Transfus Apher Sci*. 2006;35:5–17.

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